



How to change your battery

CamTAD Office
London Road Community Hospital
London Road
Derby DE1 2QY

Direct Tel: 01332 258143

E-mail: camtad@camtad.org
Website: www.camtad.org

Step one:

Your hearing aid battery lasts about a week.

Switch it off by opening the up the battery door



Step two:

Remove the battery from the hearing aid.

Batteries can be recycled at most shops and super-markets



Step three:

Remove a battery from the pack and place in the battery door.

The orange sticker should be on top.



Step four:

Remove the sticker and wait a few minutes before closing the battery door.



Supported by:

Southern Derbyshire Clinical Commissioning Group & University Hospitals of Derby & Burton NHS Foundation Trust

Patron – Professor Owen Judd, BMedSci (Hons), BM, BS, MRCPS (Glas), FRCS (Ed), FRCS (ORL-HNS), FRSA, PGDipClinEd, DCH

Charity Registration No. 517590