



How to check if a hearing aid is working

CamTAD Office

London Road Community Hospital
London Road
Derby DE1 2QY

Direct Tel: 01332 258143

E-mail: camtad@camtad.org
Website: www.camtad.org

1.) Switch the hearing aid on by closing the battery door



2.) Cup the hearing aid in both hands



3.) If the aid whistles then the aid is working. If it is not working then it may have a faulty battery or need re-tubing.

Supported by:

Southern Derbyshire Clinical Commissioning Group & University Hospitals of Derby & Burton NHS Foundation Trust

Patron – Professor Owen Judd, BMedSci (Hons), BM, BS, MRCPs (Glas), FRCS (Ed), FRCS (ORL-HNS), FRSA, PGDipClinEd, DCH

Charity Registration No. 517590